
Product Datasheet

RNA polymerase II CTD repeat YSPTSPS (phospho S2) Antibody (BF750) (orb1579240)

Description

RNA polymerase II CTD repeat

(BF750)

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

BF750

Tested

FC, IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human RNA polymerase II CTD around the phosphorylation site of Ser2

Target

 RNA polymerase II CTD repeat
YSPTSPS (phospho S2)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID
P24928
Dilution Range

Flow-Cyt=1ug/Test,IF=1:100-500

Expiration Date

12 months from date of receipt.