

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized Antibody (orb1579758)

Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid, WT & Ormicron Recognized
Preservatives	0.01M PBS(pH7.4)
Form/Appearance	Liquid
Concentration	> 1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P0DTC2
Dilution Range	WB=1:500-2000,ELISA=1:5000- 10000
Expiration Date	12 months from date of receipt.