

Product Datasheet

Phospho-Cardiac Troponin I (Ser23 + Ser24) Antibody (BF750) (orb1581921)

Description	Phospho-Cardiac Troponin I (Ser23 + Ser24) Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	BF750
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Troponin I around the phosphorylation site of Ser23/24
Target	Phospho-Cardiac Troponin I (Ser23 + Ser24)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P19429
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.