
Product Datasheet

**phospho-ATM (Ser1981) Antibody (FITC)
(orb1583792)**

Description

phospho-ATM (Ser1981) Antibody

Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ATM around the phosphorylation site of Ser1981 EG(p-S)Q
Target	ATM
Preservatives	Liquid
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	370 kDa
Uniprot ID	Q62388
Dilution Range	ICC/IF=1:50-100, IF=1:50-100, IHC-F=1:50-100, IHC-P=1:50-100, WB=1:300-1000
Expiration Date	12 months from date of receipt.