
Product Datasheet

**phospho-AKT1+AKT2+AKT3 (Tyr315+316+312)
Antibody (BF750) (orb1586379)**

Description

phospho-AKT1+AKT2+AKT3

(BF750)

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

BF750

Tested

FC, IF

Applications**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human AKT1 around the phosphorylation site of Tyr315

Target

phospho-AKT1+AKT2+AKT3 (Tyr315+316+312)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID**P31749****Dilution Range**

Flow-Cyt=1µg/Test,IF=1:100-500

Expiration Date

12 months from date of receipt.