

Product Datasheet

phospho-JAK2 (Tyr1007+Tyr1008) Antibody (BF750) (orb1586866)

Description	phospho-JAK2 (Tyr1007+Tyr1008) Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat, Zebrafish
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human JAK2 around the phosphorylation site of Tyr1007/1008
Target	phospho-JAK2 (Tyr1007+Tyr1008)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O60674
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.