

Product Datasheet

phospho-JAK2 (Tyr1007+Tyr1008) Antibody (BF680) (orb1586867)

| | |
|----------------------------|---|
| Description | phospho-JAK2 (Tyr1007+Tyr1008) Antibody (BF680) |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat, Zebrafish |
| Conjugation | BF680 |
| Tested Applications | FC, ICC, IF |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human JAK2 around the phosphorylation site of Tyr1007/1008 |
| Target | phospho-JAK2 (Tyr1007+Tyr1008) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | O60674 |
| Dilution Range | Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100 |
| Expiration Date | 12 months from date of receipt. |