

Product Datasheet

Phospho-PLCG 2 (Tyr1217) Antibody (BF680) (orb1588463)

| | |
|----------------------------|---|
| Description | Phospho-PLCG 2 (Tyr1217) Antibody (BF680) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | BF680 |
| Tested Applications | FC, IF |
| Immunogen | KLH conjugated Synthesised phosphopeptide derived from human PLC gamma 2 around the phosphorylation site of Tyr1217 |
| Target | Phospho-PLCG 2 (Tyr1217) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P16885 |
| Dilution Range | Flow-Cyt=1µg /Test,IF=1:100-500 |
| Expiration Date | 12 months from date of receipt. |