
Product Datasheet

Phospho-Estrogen Receptor alpha (Ser104 + Ser106) Antibody (BF680) (orb1590124)

Description

Phospho-Estrogen Receptor

	Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	BF680
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human ER alpha around the phosphorylation site of Ser104/106
Target	Phospho-Estrogen Receptor alpha (Ser104 + Ser106)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03372
Dilution Range	Flow-Cyt=1µg /Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.