
Product Datasheet

Phospho-Estrogen Receptor alpha (Ser104 + Ser106) Antibody (BF647) (orb1590125)

Description

Phospho-Estrogen Receptor

Antibody (BF647)

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

BF647

Tested Applications

FC, ICC, IF

Immunogen

KLH conjugated Synthesised phosphopeptide derived from human ER alpha around the phosphorylation site of Ser104/106

Target

Phospho-Estrogen Receptor alpha (Ser104 + Ser106)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID**P03372****Dilution Range**

Flow-Cyt=1µg /Test,IF=1:100-500,ICC=1:100

Expiration Date

12 months from date of receipt.