

Product Datasheet

Phospho-PTEN (Ser385) Antibody (BF750) (orb1591405)

Description	Phospho-PTEN (Ser385) Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human PTEN around the phosphorylation site of Ser385
Target	Phospho-PTEN (Ser385)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60484
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.