

Product Datasheet

phospho-beta Catenin (Tyr333) Antibody (BF680) (orb1592776)

Description	phospho-beta Catenin (Tyr333) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human beta Catenin around the phosphorylation site of Tyr333
Target	phospho-beta Catenin (Tyr333)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35222
Dilution Range	Flow-Cyt=1µg/Test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.