
Product Datasheet

Phospho-Beta-Catenin (Ser33 + Ser37) Antibody (BF750) (orb1593647)

Description

Phospho-Beta-Catenin (Ser33 +

Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Beta-Catenin around the phosphorylation site of Ser33/37
Target	Phospho-Beta-Catenin (Ser33 + Ser37)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P35222
Dilution Range	Flow-Cyt=0.2µg /test,IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.