

---

## Product Datasheet

### Phospho-Beta-Catenin (Ser33 + Ser37) Antibody (BF350) (orb1593654)

**Description**

Phospho-Beta-Catenin (Ser33 +

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	BF350
<b>Tested Applications</b>	FC, ICC, IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human Beta-Catenin around the phosphorylation site of Ser33/37
<b>Target</b>	Phospho-Beta-Catenin (Ser33 + Ser37)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P35222</a>
<b>Dilution Range</b>	Flow-Cyt=0.2µg /test,IF=1:100-500,ICC=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.