

Product Datasheet

Beta tubulin (Loading Control) Antibody (BF750) (orb1593759)

Description	Beta tubulin (Loading Control) Antibody (BF750)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	BF750
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human Beta tubulin
Target	Beta tubulin (Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P07437
Dilution Range	IF=1:200-500,ICC=1:100
Expiration Date	12 months from date of receipt.