

Product Datasheet

GPCR ORF4/GPCR GPR137 Antibody (BF680) (orb1595584)

| | |
|----------------------------|--|
| Description | GPCR ORF4/GPCR GPR137 Antibody (BF680) |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | BF680 |
| Tested Applications | ICC, IF |
| Immunogen | KLH conjugated synthetic peptide derived from human GPCR ORF4/GPCR GPR137 |
| Target | GPCR ORF4/GPCR GPR137 |
| Preservatives | Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4 |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q96N19 |
| Dilution Range | IF=1:100-500,ICC=1:100-500 |
| Expiration Date | 12 months from date of receipt. |