

Product Datasheet

Synaptotagmin 1/SYT1 Antibody (Biotin) (orb1596285)

| | |
|----------------------------|--|
| Description | Synaptotagmin 1/SYT1 Antibody (Biotin) |
| Species/Host | Rabbit |
| Reactivity | Rat |
| Conjugation | Biotin |
| Tested Applications | ELISA, ICC, IHC-Fr, IHC-P |
| Immunogen | KLH conjugated synthetic peptide derived from human Synaptotagmin 1/SYT1 |
| Target | Synaptotagmin 1/SYT1 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P21579 |
| Dilution Range | IHC-P=1:100-500,IHC-F=1:100-500,ICC=1:100-500,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |