

Product Datasheet

FMDV VP1 protein (type O) Antibody (Biotin) (orb1596327)

Description	FMDV VP1 protein (type O) Antibody (Biotin)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Recombined FMDV VP1 protein (type O)
Target	FMDV VP1 protein (type O)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P03305
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.