

Product Datasheet

H5N6-NA Antibody (AP) (orb1596779)

Description	H5N6-NA Antibody (AP)
Species/Host	Mouse
Reactivity	Virus
Conjugation	AP
Tested Applications	WB
Immunogen	Recombinant H5N6-NA protein, His
Target	H5N6-NA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	A0A240EUK5
Expiration Date	12 months from date of receipt.