

Product Datasheet

ASFV p30 Antibody (Biotin) (orb1597137)

Description	ASFV p30 Antibody (Biotin)
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombined ASFV p30 protein
Target	ASFV p30
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P34204
Expiration Date	12 months from date of receipt.