

Product Datasheet

mScarlet Antibody (HRP) (orb1597251)

Description	mScarlet Antibody (HRP)
Species/Host	Mouse
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinantm mScarlet proteins
Target	mScarlet
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9U6Y8
Expiration Date	12 months from date of receipt.