

Product Datasheet

PIK3R1 (animal-free) Antibody (HRP) (orb1597265)

Description	PIK3R1 (animal-free) Antibody (HRP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant mouse PI3K p85 Protein
Target	PIK3R1 (animal-free)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P27986
Expiration Date	12 months from date of receipt.