

Product Datasheet

PPV-VP2 Antibody (Biotin) (orb1597660)

| | |
|----------------------------|--|
| Description | PPV-VP2 Antibody (Biotin) |
| Species/Host | Mouse |
| Reactivity | Virus |
| Conjugation | Biotin |
| Tested Applications | WB |
| Immunogen | Recombinant PPV VP2 protein |
| Target | PPV-VP2 |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Form/Appearance | Liquid |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Monoclonal |
| Uniprot ID | P18546 |
| Dilution Range | WB=1:500-2000,ELISA=1:5000-10000 |
| Expiration Date | 12 months from date of receipt. |