

Product Datasheet

HPV16 E7 Antibody (BF750) (orb1597679)

Description	HPV16 E7 Antibody (BF750)
Species/Host	Rabbit
Conjugation	BF750
Tested Applications	FC, IF
Immunogen	KLH conjugated synthetic peptide derived from HPV16 E7 protein
Target	HPV16 E7
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P03129
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.