

---

## Product Datasheet

**Phospho-PDGFR $\alpha$  (Tyr849)/PDGFR $\beta$  (Tyr857)  
Antibody (BF680) (orb1598690)**

**Description**Phospho-PDGFR $\alpha$ 

Antibody (BF680)

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

BF680

**Tested Applications**

ICC, IF

**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human PDGF Receptor  $\alpha$  around the phosphorylation site of Tyr849

**Target**

Phospho-PDGFR $\alpha$  (Tyr849) / PDGFR $\beta$  (Tyr857)

**Preservatives**

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

**Form/Appearance**

Liquid

**Concentration**

1mg/ml

**Storage**

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID****P16234****Dilution Range**

IF=1:100-500,ICC=1:100

**Expiration Date**

12 months from date of receipt.