

Product Datasheet

GAPDH (Loading Control, animal-free) Antibody (Biotin) (orb1599531)

Description	GAPDH (Loading Control, animal-free) Antibody (Biotin)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Biotin
Tested Applications	IHC-P, WB
Immunogen	Full length GAPDH protein purified from rabbit muscle
Target	GAPDH (Loading Control, animal-free)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P04406
Dilution Range	WB=1:1000-10000,IHC-P=1:100-500
Expiration Date	12 months from date of receipt.