

---

## Product Datasheet

### **SARS-CoV-2 (2019-nCoV) Nucleocapsid Antibody (HRP) (orb1599830)**

**Description**

SARS-CoV-2 (2019-nCoV)

<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Recombinant SARS-CoV-2 Nucleocapsid Protei
<b>Target</b>	SARS-CoV-2 (2019-nCoV) Nucleocapsid
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	46 kDa
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Dilution Range</b>	WB=1:500-2000,ELISA=1:5000-10000
<b>Expiration Date</b>	12 months from date of receipt.