

Product Datasheet

SARS-CoV-2 (2019-nCoV) Nucleocapsid Antibody (Biotin) (orb1599831)

Description	SARS-CoV-2 (2019-nCoV) Nucleocapsid Antibody (Biotin)
Species/Host	Mouse
Reactivity	Virus
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protei
Target	SARS-CoV-2 (2019-nCoV) Nucleocapsid
Preservatives	1*PBS, 0.03% Proclin300.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	PODTC2
Dilution Range	WB=1:500-2000,ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.