

Product Datasheet

Phospho-MAP2 (Thr1613/1619) Antibody (BF647) (orb1601727)

Description	Phospho-MAP2 (Thr1613/1619) Antibody (BF647)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF647
Tested Applications	FC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MAP2 around the phosphorylation site of Thr1613 and Thr1619
Target	Phospho-MAP2 (Thr1613/1619)
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P11137
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.