

Product Datasheet

phospho-NFKB p65 (Ser311) Antibody (BF750) (orb1602289)

Description	phospho-NFKB p65 (Ser311) Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NFKB p65 around the phosphorylation site of Ser311
Target	phospho-NFKB p65 (Ser311)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q04206
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.