

Product Datasheet

S100A8 Antibody (HRP) (orb1602435)

Description	S100A8 Antibody (HRP)
Species/Host	Mouse
Reactivity	Human
Conjugation	HRP
Tested Applications	IHC-P
Immunogen	Recombinant human S100A8 protein
Target	S100A8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P05109
Expiration Date	12 months from date of receipt.