

Product Datasheet

Phospho-IKK alpha (Thr23) Antibody (BF555) (orb1602614)

Description	Phospho-IKK alpha (Thr23) Antibody (BF555)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF555
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human IKK alpha around the phosphorylation site of Thr23
Target	Phospho-IKK alpha (Thr23)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15111
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.