
Product Datasheet

**phospho-AMPK beta 1 + AMPK beta 2
(Ser182/Ser184) Antibody (BF750) (orb1603737)**

Description

phospho-AMPK beta 1 + AMPK

(BF750)

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

BF750

Tested Applications

FC, IF

Immunogen

KLH conjugated synthesised phosphopeptide derived from human AMPK beta 1 around the phosphorylation site of Ser182

Target

phospho-AMPK beta 1 + AMPK beta 2 (Ser182/Ser184)

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Liquid

Concentration

1mg/ml

Storage

Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Uniprot ID[Q9Y478](#)**Dilution Range**

Flow-Cyt=2ug/Test,IF=1:100-500

Expiration Date

12 months from date of receipt.