

Product Datasheet

PDCD4 Antibody (BF488) (orb1604619)

Description	PDCD4 Antibody (BF488)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF488
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human PDCD4
Target	PDCD4
Preservatives	Preservative: 15mM Sodium Azide, Constituents: 1% rAlbumin, 0.01M PBS, pH 7.4
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q53EL6
Expiration Date	12 months from date of receipt.