

---

## Product Datasheet

**phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199) Antibody (orb1605185)**

**Description**

phospho-PI3 kinase p85 alpha +

Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

BF750

**Tested Applications**

FC, IF

**Immunogen**

KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467

**Target**

phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199)

**Preservatives**

0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.

**Form/Appearance**

Liquid

**Concentration**

1mg/ml

**Storage**

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

**Note**

For research use only

**Isotype**

IgG

**Clonality**

Polyclonal

**Uniprot ID****P27986****Dilution Range**

Flow-Cyt=1ug/test,IF=1:100-500

**Expiration Date**

12 months from date of receipt.