

---

## Product Datasheet

**phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199) Antibody (orb1605189)**

Evolution Bioreagents.

phospho-PI3 kinase p85 alpha +

	Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	BF555
<b>Tested Applications</b>	FC, IF
<b>Immunogen</b>	KLH conjugated Synthesised phosphopeptide derived from human PI3 kinase p85 alpha around the phosphorylation site of Tyr467
<b>Target</b>	phospho-PI3 kinase p85 alpha + gamma (Tyr467 + Tyr199)
<b>Preservatives</b>	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<b>P27986</b>
<b>Dilution Range</b>	Flow-Cyt=1ug/test,IF=1:100-500
<b>Expiration Date</b>	12 months from date of receipt.