

Product Datasheet

Phospho-Mcl1 (Ser159 + Thr163) Antibody (BF750) (orb1605543)

Description	Phospho-Mcl1 (Ser159 + Thr163) Antibody (BF750)
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	BF750
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Mcl 1 around the phosphorylation site of Ser159/Thr163
Target	Phospho-Mcl1 (Ser159 + Thr163)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q07820
Dilution Range	Flow-Cyt=1ug/Test,IF=1:100-500,ICC=1:100
Expiration Date	12 months from date of receipt.