

Product Datasheet

phospho-PI 3 Kinase p110 delta (Tyr524) Antibody (BF680) (orb1606207)

Description	phospho-PI 3 Kinase p110 delta (Tyr524) Antibody (BF680)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC, IF
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PIK3CD around the phosphorylation site of Tyr524
Target	phospho-PI 3 Kinase p110 delta (Tyr524)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O00329
Dilution Range	Flow-Cyt=2ug/Test,IF=1:100-500
Expiration Date	12 months from date of receipt.