

## Product Datasheet

### TGF beta Receptor I Rabbit Polyclonal Antibody (PerCP-Cy5.5) (orb1608013)

|                            |                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------|
| <b>Description</b>         | TGF beta Receptor I Rabbit Polyclonal Antibody (PerCP-Cy5.5)                            |
| <b>Species/Host</b>        | Rabbit                                                                                  |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                       |
| <b>Conjugation</b>         | PerCP/Cy5.5                                                                             |
| <b>Tested Applications</b> | FC, IF                                                                                  |
| <b>Immunogen</b>           | KLH conjugated synthetic peptide derived from human TGF beta Receptor I (301-400/501aa) |
| <b>Target</b>              | TGFBR1                                                                                  |
| <b>Preservatives</b>       | 0.01M TBS (pH7.4) with 1% rAlbumin and 0.02% Proclin300.                                |
| <b>Concentration</b>       | 1mg/ml                                                                                  |
| <b>Storage</b>             | Store at 2-8°C. Protect from light. Avoid freezing.                                     |
| <b>Note</b>                | For research use only                                                                   |
| <b>Isotype</b>             | IgG                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                              |
| <b>Antibody Type</b>       | Primary Antibody                                                                        |
| <b>Dilution Range</b>      | Flow-Cyt=1ug/Test, IF=1:100-500                                                         |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                         |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)