

## Product Datasheet

### ATI-2341 (orb16090)

<b>Description</b>	ATI-2341, pepducin targeting the C-X-C chemokine receptor type 4 (CXCR4), is an allosteric agonist a...
<b>Note</b>	For research use only
<b>Formula</b>	$C_{104}H_{178}N_{26}O_{25}S_2$
<b>MW</b>	2256.82
<b>CAS Number</b>	[1337878-62-2]
<b>SMILES</b>	<chem>O=C(N[C@@H](CCCNC(N)=N)C(N[C@H](C(O)=O)CC(C)C)=O)[C@@H](NC([C@H](CCCCN)NC([C@H](CC(O)=O)NC([C@@]([C@H](O)C)([H])NC([C@H](CCSC)NC([C@H](CO)NC([C@H](CCCNC(N)=N)NC([C@H](CC(C)C)NC([C@H](CCCCN)NC([C@H](CCCCN)NC([C@H](CCC(N)=O)NC([C@@H](NC(CNC([C@H](CCSC)NC(CC</chem>
<b>Expiration Date</b>	12 months from date of receipt.