

Product Datasheet

phospho-IGF1R (Tyr1161) Antibody (orb1610449)

Description	phospho-IGF1R (Tyr1161) Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human IGF1R around the phosphorylation site of Tyr1161 (DI(p-Y)ET)
Target	phospho-IGF1R (Tyr1161)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	71/80 kDa
Uniprot ID	P08069
Expiration Date	12 months from date of receipt.