

## Product Datasheet

### Rat IgG / Cy5 Antibody [Out of stock] (orb1611054)

|                        |  |
|------------------------|--|
| <b>Description</b>     | Rat IgG / Cy5 Antibody   |
| <b>Preservatives</b>   | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol. |
| <b>Form/Appearance</b> | Liquid   |
| <b>Concentration</b>   | 1.0 mg/ml  |
| <b>Storage</b>         | Store at -20°C for one year. Avoid repeated freeze/thaw cycles.        |
| <b>Note</b>            | For research use only  |
| <b>Expiration Date</b> | 12 months from date of receipt.  |