

Product Datasheet

Rat IgG / PE-Cy5.5 Antibody [Out of stock] (orb1611055)

Description	Rat IgG / PE-Cy5.5 Antibody
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% glycerol.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.