

## Product Datasheet

# Mouse IgG2b Isotype Control Antibody [Out of stock] (orb1611064)

<b>Description</b>	Mouse IgG2b Isotype Control Antibody
<b>Preservatives</b>	0.01M PBS (pH7.4).
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.