

Product Datasheet

Mouse IgG1 / PerCP Isotype Control Antibody [Out of stock] (orb1611067)

Description	Mouse IgG1 / PerCP Isotype Control Antibody
Preservatives	0.01M PBS (pH7.4) with 0.5% rAlbumin and 0.02% Proclin300.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.