

## Product Datasheet

# Mouse IgG2b / PerCP Isotype Control Antibody [Out of stock] (orb1611069)

<b>Description</b>	Mouse IgG2b / PerCP Isotype Control Antibody
<b>Preservatives</b>	0.01M PBS (pH7.4) with 0.5% rAlbumin and 0.02% Proclin300.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1.0 mg/ml
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.