

Product Datasheet

RUNX1+RUNX2+RUNX3 Conjugated Antibody (orb1611814)

Description	RUNX1+RUNX2+RUNX3 Conjugated Antibody
Species/Host	Rabbit
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Recombinant
Source	Recombinant Rabbit Monoclonal Ab
Uniprot ID	Q01196
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.