

## Product Datasheet

# Collagen I alpha2 Conjugated Antibody (orb1612367)

<b>Description</b>	Collagen I alpha2 Conjugated Antibody
<b>Species/Host</b>	Goat
<b>Conjugation</b>	Biotin
<b>Storage</b>	Store at 4°C for 6 months.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Source</b>	Goat Polyclonal Ab
<b>Uniprot ID</b>	<b>P08123</b>
<b>Dilution Range</b>	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
<b>Expiration Date</b>	12 months from date of receipt.