

Product Datasheet

Mouse anti Rabbit IgG heavy chain(Fc)-HRP Conjugated Antibody (orb1612430)

Description	Mouse anti Rabbit IgG heavy chain(Fc)-HRP Conjugated Antibody
Species/Host	Mouse
Conjugation	Biotin
Storage	Store at 4°C for 6 months.
Note	For research use only
Clonality	Monoclonal
Source	Mouse Monoclonal Ab
Dilution Range	Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000
Expiration Date	12 months from date of receipt.