

## Product Datasheet

# POLR2A (Phospho-S2) Conjugated Antibody (orb1612450)

|                        |   |
|------------------------|---|
| <b>Description</b>     | POLR2A (Phospho-S2) Conjugated Antibody                       |
| <b>Species/Host</b>    | Rabbit  |
| <b>Conjugation</b>     | Biotin  |
| <b>Storage</b>         | Store at 4°C for 6 months.                                    |
| <b>Note</b>            | For research use only   |
| <b>Clonality</b>       | Recombinant   |
| <b>Source</b>          | Recombinant Rabbit Monoclonal Ab                              |
| <b>Uniprot ID</b>      | <b>P05187</b>   |
| <b>Dilution Range</b>  | Working with enzyme-conjugated streptavidin, 1: 50 - 1: 1,000 |
| <b>Expiration Date</b> | 12 months from date of receipt.                               |